

Note: the first temperature reading displayed is the most recent. If no buttons are pressed within 5 seconds the thermometer returns to standby mode. Over Temperature Warning Function If the measured body temperature is higher than the set value, the ear thermometer will save error messages and backlight prompts. (as per table on page 11)		Audible Alert Yes The Express Warranty Period will be a period of 12 months beginning on the date of purchase of the product evidenced by your dated sales receipt. You are required to provide proof of purchase as a condition of receiving Express Warranty services. You are entitled to a replacement product or repair of the product at our discretion according to the terms and conditions of this document. If your product is found to be faulty within the Express Warranty Period. This Express Warranty extends to the original purchaser only and is not transferable. Products distributed by Oricom are manufactured using new materials or new and used materials equivalent to new in performance and reliability. Spare parts may be new or equivalent to new. Spare parts are warranted to be free from defects in material or workmanship for thirty (30) days or for the remainder of the Express Warranty Period of the Oricom branded product in which they are installed, whichever is longer. During the Express Warranty Period, Oricom will, where possible repair and not replace the faulty product or part thereof. All component parts removed under this Express Warranty become the property of Oricom. In the unlikely event that your Oricom product has a recurring failure, Oricom may always, subject to the Competition Consumer Act 2010, at its discretion, elect to provide you with a replacement product of its choosing that is at least equivalent to your product in performance. No change to the conditions of this Express Warranty is valid unless it is made in writing and signed by an authorised representative of Oricom. This Express Warranty is provided by Oricom International Pty Ltd ABN 44 086 116 349, Unit 1, 4 Sovereign Place, South Windsor NSW 2756, herein after referred to as "Oricom". Oricom warrants that the product is free from defects in materials or workmanship during the Express Warranty Period. This Express Warranty does not extend to any product from which the serial number has been removed or was purchased outside of Australia. The benefits of this Express Warranty are in addition to other rights and remedies you may have under the Australian Consumer Law. Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality or the failure does not amount to a major failure. In the event of a minor failure, Oricom reserves the right to choose to repair or support@oricom.com.au. A Customer Support Team member will verify after troubleshooting with you. If your product qualifies under warranty, If so, they will give you a Product Return Authorisation number. We will then email a Return Authorisation form and a Repair Notice (if necessary), together with instructions on how to return the product for warranty service. Please note that if a Customer Support Team member advises that your product does not qualify for return, this warranty does not apply to your product. Products that are authorised to be returned to Oricom in Australia must include the following: - A completed Return Authorisation form - A copy of your Proof of Purchase (please keep your original copy) - The faulty product, including all accessories. Send the approved returns to: Oricom International Pty Ltd Locked Bag 658 South Windsor NSW 2756 Australia Please note that this Express Warranty excludes expenses incurred by you in returning any faulty product to us. You must bear any expense incurred (including postage, delivery, freight, transportation or insurance of the product) to return the faulty product to us, however, we will arrange delivery of the rep This product is manufactured in People's Republic of China under a quality system by: Guangzhou Bercom Medical Device Co., Ltd. No. 38, Huanchun Xi Road, Dagang Town, Nansha, 511470, Guangzhou City, Guangdong Province, P.R. CHINA IMPORTED & DISTRIBUTED IN AUSTRALIA by the sponsor- Oricom International Pty Ltd, Unit 1, 4 Sovereign Place, South Windsor 2756 AUSTRALIA.																																			
Technical Specifications <table> <tr> <td>Operational principle</td><td>Infrared radiation detection</td></tr> <tr> <td>Mode of operation</td><td>Adjustable mode</td></tr> <tr> <td>Measurement reference body site</td><td>Oral</td></tr> <tr> <td>Patient population</td><td>People of all ages</td></tr> <tr> <td>Rx/OTC Classification</td><td>OTC (Over The Counter)</td></tr> <tr> <td>Measurement site</td><td>Ear</td></tr> <tr> <td>Measuring range</td><td>34.0°C-43.0°C (93.2°F-109.4°F)</td></tr> <tr> <td>Accuracy</td><td>34.0°C-34.9°C (93.2°F-94.8°F): ±0.3°C (0.6°F) 35.0°C-42.0°C (95.0°F-107.6°F): ±0.2°C (0.4°F) 42.1°C-43.0°C (107.8°F-109.4°F): ±0.3°C (0.6°F)</td></tr> <tr> <td>Clinical repeatability</td><td>±0.2°C (0.4°F)</td></tr> <tr> <td>Batteries</td><td>DC 3V (3pcs AA dry batteries)</td></tr> <tr> <td>Measurement time</td><td>Approx. 3.5 seconds</td></tr> <tr> <td>Display mode</td><td>LCD screen</td></tr> <tr> <td>Display resolution</td><td>0.1°C (0.1°F)</td></tr> <tr> <td>Automatic power-off</td><td>Approx. 60 to 10 seconds</td></tr> <tr> <td>Memory</td><td>32 sets</td></tr> </table>		Operational principle	Infrared radiation detection	Mode of operation	Adjustable mode	Measurement reference body site	Oral	Patient population	People of all ages	Rx/OTC Classification	OTC (Over The Counter)	Measurement site	Ear	Measuring range	34.0°C-43.0°C (93.2°F-109.4°F)	Accuracy	34.0°C-34.9°C (93.2°F-94.8°F): ±0.3°C (0.6°F) 35.0°C-42.0°C (95.0°F-107.6°F): ±0.2°C (0.4°F) 42.1°C-43.0°C (107.8°F-109.4°F): ±0.3°C (0.6°F)	Clinical repeatability	±0.2°C (0.4°F)	Batteries	DC 3V (3pcs AA dry batteries)	Measurement time	Approx. 3.5 seconds	Display mode	LCD screen	Display resolution	0.1°C (0.1°F)	Automatic power-off	Approx. 60 to 10 seconds	Memory	32 sets	Batteries - The thermometer uses 2 x AA Alkaline batteries - When the battery symbol  appears on the screen and flashes, it indicates that the battery is low and needs to be replaced - If the ambient temperature is low, it will shorten the battery life. - Use alkaline batteries only - When not in use for a long time, it is recommended to remove the batteries to avoid battery leakage and damage to the product. - Note: Battery leakage not covered by our Warranty. - Disposing of batteries - To protect environment, dispose of batteries at appropriate collection sites according to national or local regulations.  Notes on handling batteries - If your skin or eyes come into contact with battery fluid, rinse the affected area with water and seek medical assistance. - Choking hazard! - Children may swallow and choke on batteries. Store batteries out of the reach of children. - If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth. - Protect batteries from excessive heat. IP22 TEMPERATURE Due to the IP22 rating, the measuring temperature above 43.0°C (109.4°F) and displays "Lo" when measuring temperature below 34.0°C (93.2°F). - Use identical or equivalent battery types only - Always replace all batteries at the same time - Do not use rechargeable batteries - Do not disassemble, split or crush the batteries.					
Operational principle	Infrared radiation detection																																				
Mode of operation	Adjustable mode																																				
Measurement reference body site	Oral																																				
Patient population	People of all ages																																				
Rx/OTC Classification	OTC (Over The Counter)																																				
Measurement site	Ear																																				
Measuring range	34.0°C-43.0°C (93.2°F-109.4°F)																																				
Accuracy	34.0°C-34.9°C (93.2°F-94.8°F): ±0.3°C (0.6°F) 35.0°C-42.0°C (95.0°F-107.6°F): ±0.2°C (0.4°F) 42.1°C-43.0°C (107.8°F-109.4°F): ±0.3°C (0.6°F)																																				
Clinical repeatability	±0.2°C (0.4°F)																																				
Batteries	DC 3V (3pcs AA dry batteries)																																				
Measurement time	Approx. 3.5 seconds																																				
Display mode	LCD screen																																				
Display resolution	0.1°C (0.1°F)																																				
Automatic power-off	Approx. 60 to 10 seconds																																				
Memory	32 sets																																				
Pack Contents 1 x Ear thermometer 1 x User manual 2 x AA Alkaline batteries 20 x Probe Covers		Troubleshooting If you have problems while using your thermometer, please refer to this guide to help resolve the problem. If the problem persists, please contact Oricom Customer service. <table> <tr> <th>Screen message</th><th>Description</th><th>Possible causes</th><th>Solutions</th></tr> <tr> <td>No display</td><td></td><td>1. System error. 2. The battery is completely depleted.</td><td>1. Turn off and then restart the thermometer. 2. Remove the battery and insert a new one. 3. Please contact Oricom Customer Support</td></tr> <tr> <td></td><td>Battery symbol flashing, can't measure temperature</td><td>The current is too low to perform temperature measurement</td><td>Replace batteries.</td></tr> <tr> <td></td><td>Battery symbol flashing, last battery level flashing</td><td>Batteries are low, but still operate correctly</td><td>Suggest replacing batteries as soon as possible.</td></tr> <tr> <td></td><td>Probe cover symbol flashing, animated style</td><td>1. No probe cover installed. 2. The probe cover is not properly installed.</td><td>1. Install probe cover. 2. Remove the probe cover and reinstall.</td></tr> <tr> <td>Display "Er1"</td><td>Ambient temperature is not within the measuring range (10°C to 40°C)</td><td>Use product within the ambient temperature operating range (10°C to 40°C)</td><td></td></tr> <tr> <td>Display "Hi"</td><td>The measured temperature exceeds the measurement range (34.0°C to 43.0°C)</td><td>1. Check if the probe lens is dirty. Please clean the probe lens and operate again. 2. Please ask the subject to maintain a stable posture and measure again according to the correct operating method. 3. Place the product at normal temperature (10°C to 40°C) and wait for 30 minutes.</td><td></td></tr> </table>		Screen message	Description	Possible causes	Solutions	No display		1. System error. 2. The battery is completely depleted.	1. Turn off and then restart the thermometer. 2. Remove the battery and insert a new one. 3. Please contact Oricom Customer Support		Battery symbol flashing, can't measure temperature	The current is too low to perform temperature measurement	Replace batteries.		Battery symbol flashing, last battery level flashing	Batteries are low, but still operate correctly	Suggest replacing batteries as soon as possible.		Probe cover symbol flashing, animated style	1. No probe cover installed. 2. The probe cover is not properly installed.	1. Install probe cover. 2. Remove the probe cover and reinstall.	Display "Er1"	Ambient temperature is not within the measuring range (10°C to 40°C)	Use product within the ambient temperature operating range (10°C to 40°C)		Display "Hi"	The measured temperature exceeds the measurement range (34.0°C to 43.0°C)	1. Check if the probe lens is dirty. Please clean the probe lens and operate again. 2. Please ask the subject to maintain a stable posture and measure again according to the correct operating method. 3. Place the product at normal temperature (10°C to 40°C) and wait for 30 minutes.							
Screen message	Description	Possible causes	Solutions																																		
No display		1. System error. 2. The battery is completely depleted.	1. Turn off and then restart the thermometer. 2. Remove the battery and insert a new one. 3. Please contact Oricom Customer Support																																		
	Battery symbol flashing, can't measure temperature	The current is too low to perform temperature measurement	Replace batteries.																																		
	Battery symbol flashing, last battery level flashing	Batteries are low, but still operate correctly	Suggest replacing batteries as soon as possible.																																		
	Probe cover symbol flashing, animated style	1. No probe cover installed. 2. The probe cover is not properly installed.	1. Install probe cover. 2. Remove the probe cover and reinstall.																																		
Display "Er1"	Ambient temperature is not within the measuring range (10°C to 40°C)	Use product within the ambient temperature operating range (10°C to 40°C)																																			
Display "Hi"	The measured temperature exceeds the measurement range (34.0°C to 43.0°C)	1. Check if the probe lens is dirty. Please clean the probe lens and operate again. 2. Please ask the subject to maintain a stable posture and measure again according to the correct operating method. 3. Place the product at normal temperature (10°C to 40°C) and wait for 30 minutes.																																			
Signs And Symbols  Trade mark  Indicates the need for the user to consult the instructions for use for critical cautionary information such as warnings and precautions that can, for a variety of reasons, be presented on the medical device itself.  Please read the instructions carefully.  Degree of protection against harmful ingress of water and particulate matter. 1. Protection against solid foreign objects of 12.5mm Ø and greater. 2. Protection against vertically falling water drops when EN60529, classed up to IP22.  DISPOSAL: Do not dispose of this product as unsorted municipal waste. Collection of such waste separately for special treatment is necessary.		EMC Declaration The ME EQUIPMENT or ME SYSTEM is suitable for home healthcare and clinic environments. Statement: The ear thermometer model IET800 was tested according to the recommendations of Technical Report IEC TR 60601-4-2: Medical electrical equipment – Part 4-2: Guidance and interpretation – Electromagnetic immunity: performance of medical electrical equipment and medical electrical systems. Warning: Don't use near active HF surgical equipment and the RF shielded room or an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high. Warning: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment be about to be observed to verify that they are operating normally. Warning: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation. Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the equipment IET800 including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result. All necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE with regard to electromagnetic disturbances for the accepted service life. Portable and mobile RF communication equipment may affect the performance of the IET800 avoiding strong electromagnetic interference when used, such as near mobile phones, microwave ovens, etc. Guidance and manufacturer's declaration – electromagnetic emissions and Immunity Guidance and manufacturer's declaration – electromagnetic emissions and Immunity																																			
Table 1 Guidance and manufacturer's declaration – electromagnetic emissions <table> <tr> <th>Emissions test</th><th>Compliance</th></tr> <tr> <td>RF emissions CISPR 11</td><td>Group 1</td></tr> <tr> <td>RF emissions CISPR 11</td><td>Group B</td></tr> <tr> <td>Harmonic emissions IEC 61000-3-2</td><td>Not applicable</td></tr> <tr> <td>Voltage fluctuations/flicker emissions IEC 61000-3-3</td><td>Not applicable</td></tr> </table>		Emissions test	Compliance	RF emissions CISPR 11	Group 1	RF emissions CISPR 11	Group B	Harmonic emissions IEC 61000-3-2	Not applicable	Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable	Table 3 Guidance and manufacturer's declaration – electromagnetic immunity <table> <tr> <th>Test</th><th>Test frequency (MHz)</th><th>Band</th><th>Service</th><th>Modulation</th><th>Modulation rate (Hz)</th><th>Distance (m)</th><th>EN60601-1-2 TEST LEVEL (V/m)</th></tr> <tr> <td>Radio</td><td>100 kHz to 100 MHz</td><td>380-380</td><td>TETRA 460</td><td>Pulse modulation 1.4 kHz tone</td><td>1.8</td><td>0.3</td><td>27</td></tr> <tr> <td>Radio</td><td>100 kHz to 100 MHz</td><td>380-380</td><td>GSM 1800</td><td>FM modulation 1.4 kHz tone</td><td>2</td><td>0.3</td><td>28</td></tr> <tr> <td>Radio</</td></tr></table>	Test	Test frequency (MHz)	Band	Service	Modulation	Modulation rate (Hz)	Distance (m)	EN60601-1-2 TEST LEVEL (V/m)	Radio	100 kHz to 100 MHz	380-380	TETRA 460	Pulse modulation 1.4 kHz tone	1.8	0.3	27	Radio	100 kHz to 100 MHz	380-380	GSM 1800	FM modulation 1.4 kHz tone	2	0.3	28	Radio</
Emissions test	Compliance																																				
RF emissions CISPR 11	Group 1																																				
RF emissions CISPR 11	Group B																																				
Harmonic emissions IEC 61000-3-2	Not applicable																																				
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable																																				
Test	Test frequency (MHz)	Band	Service	Modulation	Modulation rate (Hz)	Distance (m)	EN60601-1-2 TEST LEVEL (V/m)																														
Radio	100 kHz to 100 MHz	380-380	TETRA 460	Pulse modulation 1.4 kHz tone	1.8	0.3	27																														
Radio	100 kHz to 100 MHz	380-380	GSM 1800	FM modulation 1.4 kHz tone	2	0.3	28																														
Radio</																																					